

# Precise measurement of product temperature in freeze drying allows development of a scientifically justified lyo-cycle adaptable to scale and specific freeze dryer equipment features and the control of the lyo-cycle in commercial production

**Anton Mangold** IQ Mobil solutions GmbH, Industriestrasse 7, 83607 Holzkirchen, Germany

Presented by Kazuo Hayashi; Nippon Science Core & Morio Tojo; East Core Co.


## INTRODUCTION:

The precise measurement of product temperature ( $T_p$ ) is the most important parameter in rational development of lyo-cycles. The use of  $T_p$  as parameter in laboratory/pilot scale up to transfer into commercial production of different lyo-formulations will be presented. Also  $T_p$  can be applied to control the lyo-cycle of commercial production by integration of the  $T_p$  signal provided by TEMPRIS® into the PLC (SCADA - supervisory control and data acquisition) of the lyophilizer.

The application of precise measurement of  $T_p$  using a PAT tool in freeze dryers allows to scientifically justify and rationalize lyo-cycle development during all stages of lyo-cycle development from lab to commercial production leading to significantly increased process knowledge and more reliable and robust processes.  $T_p$  may be used to control the lyo-cycle by integration of the data provided by TEMPRIS into the PLC and by defining acceptance criteria for the control.


## MATERIALS & METHODS:

### Functional principle



Quartz based sensor, operating on the principle of temperature dependent resonance: after excitation by a modulated microwave signal (2.4 GHz, worldwide available ISM-band) the sensor keeps on oscillating in a temperature dependent frequency. Overlaying the sensor response with the carrier signal leads to a frequency shift from which the product temperature  $T_p$  can be derived.

### TEMPRIS® modular system for cycle development, scale up, transfer



## RESULTS & DISCUSSION (I):

### Precise measurement of $T_p$





Figure 3-83: Evaluation of edge effects using TEMPRIS sensors in a 50 mg/mL sucrose run

Tempris edge sensors are more accurate than thermocouple  $T_p$  is lower. <sup>1)</sup>


## RESULTS & DISCUSSION (II):

### Testing in dual-chambered syringe by Pfizer <sup>2)</sup>



Slide courtesy of Brian Wilbur's presentation, Process Development of a Dual-Chambered Syringe, presented at CHI PepTalk 2012.

### Transfer: „hot and cold spot“ detection <sup>3)</sup>



Lyo-cycle production scale 30 sqm Lyo, 18 shelves,

## CONCLUSION:

- As TEMPRIS® is a modular – mobile system it may be applied in all stages of lyo-cycle development from lab to commercial production.
- TEMPRIS® allows for the most precise measurement of  $T_p$  of all currently available PAT tools usable also in commercial scale aseptic lyophilization.
- Take the possibility and use data about  $T_p$  gained by TEMPRIS® as parameter to actually control the freeze drying process by placing TEMPRIS® into vials, located at the worst case positions of the lyo and by using their  $T_p$  signal as the criterion to start the next process step in the lyo-cycle.

## REFERENCES:

- 1) cf. Schneid, S. PhD Thesis 2009; p. 163
- 2) Brian Wilbur, Pfizer; CHI PepTalk ;2012
- 3) Dr. Andrea Weiland-Waibel, Explicat Pharma GmbH; 2010